



Microscopy and Cytometry Facility

Sample Preparation Guidelines (Flow Cytometry Analysis and Cell Sorting)

Sample

- Prepare single cell suspension.
- Re-suspend cells in filtered PBS buffer ($\text{Ca}^{2+}/\text{Mg}^{2+}$ free) with 0.1-1% BSA + 2.5 mM EDTA) but, if necessary, in any other medium without phenol red and serum (FCS, FBS).
- Remove cells and tissue aggregates by filtering the sample with 30-40 μm mesh filters (cell strainer). MCF can provide cell strainers.
- The sample should have particles density of 5-10 million/ml. If less cells are available, put then in a minimum volume of 500 μl .
- Bring the samples in any leak-proof tubes or in a 5ml round-bottom polypropylene tube. MCF can provide 5ml round-bottom tubes.
- Provide additional buffer in case the sample for sorting needs to be diluted.
- Provide collection tubes/plates one-fifth filled with collection buffer or medium. Note that the collected sample will be diluted with PBS from the instrument sheath fluid.

Tips

- To improve cells viability, all preparation steps should be performed on ice unless otherwise specified in your protocol.
- Use minimal speed to sediment cells (start with 300xg). Avoid vigorous vortexing, drying of the cell pellet and introduction of air bubbles during sample preparation.
- To exclude dead cells from the analysis or sorting it is HIGHLY recommended to use the so called live/dead dyes (DNA-binding or amine-reactive dyes) as markers of dead cells.
- DNA released from dead cells during cells dissociation causes stickiness. To reduce cells aggregation in your samples, incubate the cells in the presence of DNase I (100 $\mu\text{g}/\text{ml}$ with 5 mM MgCl_2).

Gating and Compensation controls

- Each single-color or multi-color experiment requires appropriate staining controls to properly set the population gates. Please provide:
 - Negative control - unstained sample (to detect autofluorescent background),
 - Positive controls - single stained samples (to measure positive signals and eliminate spillover to other detectors),
 - FMO controls (Fluorescence Minus One) for multi-color experiments,